

Optical Scattering Processes for Tissue Diagnosis


Elastic Scattering Processes: Reflectance, Structural OCT


Inelastic Scattering Processes: Fluorescence, Spontaneous Raman

Stimulated Coherent Inelastic Scattering Processes:
 Stimulated emission, Stimulated Raman scattering
 Four-wave mixing processes, SHG, THG
 Coherent Anti-Stokes Raman Scattering (CARS)


• Coherent Anti-Stokes Raman Scattering (CARS) produces
 measurable signal with sufficient power


• Take advantage of the coherent nature of CARS and use
 interferometry to provide both amplitude and phase information
 for molecular diagnostics








Nonlinear Interferometric Vibrational Imaging (NIVI) Interferometric detection based on a coherent reference signal Amplitude and Phase determination (Im[χ(3)] = Raman spectra) NIVI is linear in molecular concentration Raman lineshape is free of non-resonant background Enables quantitative spectral decomposition for analyzing mixtures or differentiating subtle spectral features between molecules Significantly faster than Raman spectroscopy, and enables imaging









































































































Driving Hypothesis 1. Sparse cancer cells release extracellular vesicles to modulate their microenvironment inside the tissue. Numerous normal cells of the tissue accept the vesicles to switch their metabolism from energy production to biomass production.


- 3. The switched metabolism induces a structurally modified microenvironment (tumor organ) that promote the proliferation, invasion, and metastasis of the cancer cells.

















































Conclusions

- o Broadband coherent fiber light source
- o Optimal pulse shaping
- o Label-free multi-modal multiphoton imaging and molecular profiling
- o Multiphoton cancer biomarkers
- Optical label-free molecular imaging for cancer diagnosis and profiling the tumor microenvironment
- o Label-free in situ molecular histopathology and surgical guidance

Liu et al. Label-free molecular profiling for identification of biomarkers in mor tumor development using multimodal multiphoton imaging. Submitted 2015





Future Directions

- o Fast-scanning, image in vivo dynamics
- o Developing mathematical formulism to quantify inter-related multimodal changes
- Isolate and characterize microvesicles
- o Identify other recurring biomarkers
- o Correlate pre-clinical with clinical findings
- o Engineer portable system
- o Intraoperative imaging during breast cancer surgery
- o Investigate pharmacological interventions





